

Feng Wang


On behalf of the HGCal IHEP group

27/11/2021



Structure of CMS endcap calorimeter









CMS HGCal module





- Absorber: PCB (FR4) baseplate for CE-H & CuW baseplate for CE-E
- Sensor: 8 inch with120µm, 200µm and 300µm depletion depth for different radiation level
- PCB: Hexaboard with front-end electronics ASIC (HGCROC)
- Assemble module by glue in each layer




- Wire bonding for signal read out, guard ring and BV
- Protect wire by encapsulation













• Module assembly lab

- > Clean room with temperature: $21 \pm 1^{\circ}$ C and humidity: $45\% \pm 10\%$
- ≻ Key equipment: Gantry, bonder, OGP, pull tester, mini gantry, teststand, microscope, etc.

• Selection of components

- > Flatness of baseplate and sensor: less than $\pm 50 \mu m$
- \succ No fatal scratches found on the sensor
- > Acceptable noise level for front-end board (Hexaboard)
- Key indicators for QC in the module assembly process
 - \succ Glue coverage area greater than 70% for heat dissipation
 - > Thickness of epoxy glue layer in the assembly procedure: $100\mu m \pm 25\mu m$
 - > Edge-to-edge alignment less than: $\pm 30 \mu m$
 - > Pulling force of a single bonding wire >8g
 - > Single module (LD) with 192 channels requires more than 600 wires
 - > All bonding wires are encapsulated
- Module testing

Good noise and IV performance before and after assembly into modules


We have met the above requirements at IHEP MAC



IHEP HGCal lab



- 140m² dedicated clean room (Class 1000) at the ready:
 - Continuous operation for one year
 - Temperature and humidity have good performance
 - Compressed/vacuum/dry air pipe installed
 - Key equipment is ready
 - Training over





Dispense

Wire bonding

Pull test & encapsulation


Module test



HGCal module automatic assembly procedure



- Check components before assembling
 - ➤ Thickness, size, etc. measurement for baseplate by OGP
 - Scratch inspection on sensor by OGP
 - Pedestal, noise, etc. testing for hexaboard by testand
- Single module automatic assembly
 - Step 1: the baseplate, sensor and hexaboard assembly on the gantry
 - Step 2: wire bonding for read-out, guard ring and BV are operated on the bonder
 - Step 3: test module with wire bonding and without
 - Step 4: encapsulation for protecting bonding wire is carried out on the mini gantry
 - Step 5: teststand is for the full tests of the noise and IV performance of module (after assembly)






High precision automatic assembly




7

- Automated assembly on the gantry
 - Calibrate the components position by CCD
 - Dispense on the component
 - Pick up the specific tooling by vacuum
 - Pick and place sensor (hexaboard) on baseplate (sensor)
 - Curing for 24 hours
- Performance
 - ➢ Glue thickness meet the requirement
 - ➢ Glue cover area lager than 70%




Curing



Assembly Hexaboard on sensor

Glue thickness: 1.369-0.300 - 0.965= 0.104mm







High speed automatic wire bonding




- Wire bonding
 - > 4 wires (23 μ m) can be bonded per second
 - Specific tooling holds module on the bonder
 - \succ Front side bonding is for read-out and guard ring
 - Back side bonding is for bias voltage
 - > 3 wires for each pixel
- Performance
 - > Pulling forces of a single wire greater than 8g





Read out bonding




Guard ring bonding







- ➤ Glue: Sylgard 186 1:10 (mass) which is a clear silicon glue with very high radiation tolerance
- Centrifuge is used to remove air bubbles
- > Automatic glue dispense for encapsulation on the mini gantry








First 8 inch HGCal Module @IHEP




- Test leakage, Pedestal, Noise, etc.
- First test beam for 8 inch HGCal module at CERN


















First 8 inch HGCal module assembly success at IEHP

- > All assembly procedure has been video recorded step by step and shared with CMS HGCal group
- > The thickness control, wire pulling force and encapsulation of this module mainly meet the requirements
- > The electronics test of this module is normal, and detects the signal in the beam test at CERN
- > We have mainly mastered the key technology of the single HGCal module assembly
- ➤ Experts from CMS HGCal group express their congratulations to IHEP



- We have done an other 8 inch HGCal module at IHEP
 - \succ The performance is similar with the first one
 - \blacktriangleright We sent this module to CERN and tested the beam test on Oct. 2nd, 2021
- **Prospects**
 - Optimize the module assembly process
 - Product assembly tooling to match the mass production





顾煜栋

华慧玲

王储

张华桥

刘勇

王峰

陶军全

程华杰

王泽炳

Thankful to close collaborations with HGCAL teams, esp. CERN, UCSB and NT

iogo

100

郭佳林

余涛哲

苑超辰