

High-granularity Crystal ECAL for CEPC

Fangyi Guo on behalf of the CEPC Calorimeter Working Group IHEP CAS

CEPC International Workshop, Nov. 2021



中国科学院高能物理研究所

Institute of High Energy Physics Chinese Academy of Sciences

Outline

Introduction.

- PFA performance with crystal cube ECAL:
 - Separation power power with 2 incident particles.
 - Higgs benchmark @ CEPC $H \rightarrow \gamma \gamma$ (2 photons), $H \rightarrow gg$ (2 jets).
- A new proto-PFA software for crystal bar ECAL.
- Technical developments: crystals and SiPMs:
 - Uniformity: Geant4 simulation vs. Lab measurement.
 - Crystal measurement: energy resolution.



Introduction

- Future high energy lepton collider (e.g. CEPC):
 - Excellent high-energy jet energy measurement (3~4% resolution @100GeV) for Higgs and EW study.
 - Precise γ/π^0 separation for flavor physics and BSM.
- Particle-flow Approach (PFA):
 - Measure jet by its components: 60% charged particles, 30% photons, 10% neutral hadrons.

• Final resolution:
$$\sigma_{Jet} = \sqrt{\sigma_{track}^2 + \sigma_{EM}^2 + \sigma_{Had}^2 + \sigma_{confusion}^2}$$

Component	Detector	Energy Fraction	Energy Resolution	Jet Energy Resolution
Charged Particles (X^{\pm})	Tracker	$\sim 0.6 E_J$	_	—
Photons (γ)	ECAL	$\sim 0.3 E_J$	$0.15\sqrt{E_{\gamma}}$	$0.08\sqrt{E_J}$
			$0.03 \sqrt{E_{\gamma}}$	$0.016\sqrt{E_J}$
Neutral Hadrons (h^0)	HCAL	$\sim 0.1 E_J$	$0.55 \sqrt{E_{h^0}}$	$0.17 \sqrt{E_J}$

Jet E res.	W/Z sep	
perfect	3.1 σ	
2%	2.9 σ	
3%	2.6 σ	
4%	2.3 σ	
5%	2.0 σ	
10%	1.1 σ	



Introduction

- PFA requirement: Hardware + Software
 - Distinguish showers in calorimeter ⇒ high granularity ECAL/HCAL.
 - Minimize transverse spread of EM shower \Rightarrow small Moliere radius R_M
 - \Rightarrow SiW sampling ECAL in ILD.
 - Separate EM and Hadronic showers longitudinally \Rightarrow large λ_I/X_0 ratio.
- Crystal ECAL:
 - Homogeneous structure \Rightarrow energy resolution $\sim 3\%/\sqrt{E} \oplus 1\%$.
 - Capability to trigger single photons \Rightarrow precise γ/π^0 reconstruction.
 - High sensitivity to low energy particles
 essential for flavor physics.
 - Larger Moliere radius
 larger probability of shower overlap.
 - Smaller $\lambda_I / X_0 \Rightarrow$ larger probability of hadronic shower in ECAL.

Material	X_0 /cm	R_M /cm	λ_I /cm	λ_I/X_0
W	0.35	0.93	9.6	27.4
BGO	1.12	2.23	22.8	20.3
Ratio	3.2	2.4	2.4	0.74







2 major designs for crystal ECAL

Design 1: high-granularity crystal cubes



- Fine segmentation of crystal cube $(1 \times 1 \times 1 cm^3)$, single-ended readout with SiPM.
- Compatible with PFA.
- Focus on PFA performance studies with ArborPFA.

Design 2: cross arranged crystal bars



- Long crystal bars: 1×1×40 cm³, doubleside readout with SiPM.
- Crossed arrangement in adjacent layers + timing at 2 sides for positioning.
- Super cell module: 40×40 cm²
- Save #channels and minimize dead materials.
- Focus on reconstruction algorithm
 development
 2021/11/9



- Confusion term in jet energy resolution $\sigma_{confusion}$:
 - Pattern recognition: distinguish nearby clusters from different particles.
 - But how often will this confusion happen in detector?
- Physics topology for 4-jet event in CEPC: $ee \rightarrow ZH \rightarrow qqgg$.
 - Multiplicity in a $40 \times 40 \ cm^2$ tower: ≤ 2 particles for 90% cases.
 - Separation between $\gamma + \gamma$ and $\gamma + \pi^{\pm}$ can be a mark of σ_{conf} .



104 visible final state particles, 260 pairs with distance < 400mm Hottest cell: 10 hits, E / 1/s = 32.8GeV / 240.0GeV = 13.7% Particle hit position



Outline

- Introduction.
- PFA performance with crystal cube ECAL:
 - Separation power power with 2 incident particles.
 - Higgs benchmark @ CEPC $H \rightarrow \gamma \gamma$ (2 photons), $H \rightarrow gg$ (2 jets).
- A new proto-PFA software for crystal bar ECAL.
- Technical developments: crystals and SiPMs:
 - Uniformity: Geant4 simulation vs. Lab measurement.
 - Crystal measurement: energy resolution.



- Diphoton separation in ArborPFA: overlap between two EM shower.
 - Dominant by Moliere radius of BGO: 2.23cm.
 - Solution: first use a higher energy threshold to extract EM shower core for separation, then use low threshold to recover the energy (WIP).
 - With a proper energy threshold crystal ECAL can reach similar separation power with SiW ECAL. gamma/gamma Separation Efficiency



gamma/gamma Separation Efficiency

200

250

Distance / mm



Particle gun event of 2 5GeV photons. Scan the angle(distance) between two photons, check the successful reconstruction efficiency.



- $\gamma + \pi^+$ separation in ArborPFA: much more complex.
 - BGO vs. SiW: larger transverse spread (R_M) and longitudinal overlap probability (λ_I/X_0) .
 - High threshold may lose hadronic shower fractions.
 - Parameter tuning in ArborPFA: cluster merging.





Particle gun event of 5GeV γ + 10GeV π^+ . Scan the angle(distance) between γ and π^+ , check the successful reconstruction efficiency.



- Physics benchmark: full simulation of ZH process:
 - $e^+e^- \rightarrow ZH \rightarrow \nu\nu\gamma\gamma$ process @ 240GeV: BMR=1.1%
 - Compensation and correction in ECAL barrel gaps has not been finished in crystal ECAL, hits near gaps are excluded.
 - $e^+e^- \rightarrow ZH \rightarrow \nu\nu gg$ di-jet process @ 240GeV:
 - Tuned the parameters in ArborPFA: digitization threshold, bushconnect parameters, etc.
 - Get significant improvement from 4.5% @<u>Yangzhou</u> to present 4.0%, get closer to SiW result 3.8%.





Outline

- Introduction.
- PFA performance with crystal cube ECAL:
 - Separation power power with 2 incident particles.
 - Higgs benchmark @ CEPC $H \rightarrow \gamma \gamma$ (2 photons), $H \rightarrow gg$ (2 jets).
- A new proto-PFA software for crystal bar ECAL.
- Technical developments: crystals and SiPMs:
 - Uniformity: Geant4 simulation vs. Lab measurement.
 - Crystal measurement: energy resolution.



Proto-PFA for crystal bar ECAL

- Crystal Bar ECAL design: more tracker-like structure.
 - Less #channels and materials in ECAL.
 - Ambiguity problem from 2D measurement (ghost hit)
 - Identification of energy deposits from particles (confusion)
- A special software for this ECAL design:
 - basic unit: double layer with crossed bars.
 - Combine 2 layers to mimic a high-granularity ECAL.
 - Use time + similar energy in adjacent layer for crosslocation.





Proto-PFA for crystal bar ECAL

- Reconstruction flow:
 - 1D: clustering and energy splitting.
 - Cluster the neighbor fired bars.
 - Use EM shower profile to split two nearby showers.
 - 2D: Energy + time matching in 2 adjacent layers for ghost hit removal
 - Define a χ^2 with both energy and timing info.
 - Reject the wrong combination with χ^2 .
 - 3D: Cluster ID + cone clustering.
 - Use lateral moment $LAT = \frac{\sum_{i=3}^{N} E_i r_i^2}{\sum_{i=3}^{N} E_i r_i^2 + E_1 r_0^2 + E_2 r_0^2}$ for cluster ID
 - longitudinally cone clustering for MIP/EM/Hadronic showers.
 - WIP: cluster merging, track-cluster matching, etc.



Proto-PFA for crystal bar ECAL

- Separation: very preliminary result.
 - Diphoton separation: two 5GeV photons



• $\gamma + \pi^{-}$ separation: 5GeV $\gamma + 10$ GeV π^{-} . WIP.







Outline

- Introduction.
- PFA performance with crystal cube ECAL:
 - Separation power power with 2 incident particles.
 - Higgs benchmark @ CEPC $H \rightarrow \gamma \gamma$ (2 photons), $H \rightarrow gg$ (2 jets).
- A new proto-PFA software for crystal bar ECAL.
- Technical developments: crystals and SiPMs:
 - Uniformity: Geant4 simulation vs. Lab measurement.
 - Crystal measurement: energy resolution.



Crystal uniformity measurement

4775

33.36

14.3 0

7.725 / 10

 207.7 ± 6.2

 41.29 ± 0.24

 6.253 ± 0.430

80

0.6556

0

- BGO response simulation with Geant4 10.7:
 - 400mm BGO crystal bar, transverse 1cm²
 - 662keV gamma from Cs-137
 - Varying Cs-137 position.

240 f

220

200

180

160

140

120F

100

80

60

40

20 F

JYJ_

10

20

30

Number

fit 662keV photon peak to get #photons

662keV photopeak

40

HistDetectedPhoton



• Generally good response uniformity expected in G4 simulation.

70

Detected Photon

Entries

Mean

Std Dev

Underflow

Overflow

Constant

 χ^2 / ndf

Prob

Mean

Sigma

60

(gamma hitting the center)

50

Crystal uniformity measurement

- BGO response in lab:
 - Setup: 400mm long BGO crystal (with ESR foil) and Cs-137 source.
 - The same configuration as the simulation.
 - Trends are not significant enough due to the systematic difference between 2 SiPMs.
 - Work plan: to use optical grease to improve the crystal-SiPM coupling and reproducibility







Crystal response in Lab

Photon energy resolution with impacts of:

Surfaces: polished/ground.





100000

553.6

248.3

0.381

99871

534.6

210.3

5034

20.11/30

 1648 ± 10.4

 662.9 ± 1.1

86.77 ± 1.19

1200

Energy [keV]

0.9137

298

182

1.128e+04

49.31 / 47

 1314 ± 8.1

 661.6 ± 0.9

 124.7 ± 1.4

Energy [keV]





2021/11/9

Energy Resolution (E.R.) = $2.355 \times \frac{\sigma}{mean}$, defined as FWHM



- PMT has better acceptance (full coverage of crystal transverse area) than SiPM, to be updated with larger SiPMs
- Further comparisons will be done with simulation

Summary

- Steady progress to address key issues.
- Performance studies with crystals using Arbor-PFA:
 - Separation power of close-by particles.
 - Performance studies with Higgs events: closer to SiW.
- Developing a new proto-PFA software for crystal ECAL:
 - Traditional PFA: fine granularity + small R_M + less hits (sampling) for separation.
 - Crystal PFA: precise energy (homogeneous) + shower profile for separation.
 - Key issue: ghost hit and confusion problem can be solved. Preliminary result is promising.
 - Many details still need optimization.
- Technical developments:
 - Good uniformity with long crystal bar, and experiment can match Geant4 simulation.
 - BGO crystal shows better than 20% energy resolution to 662keV photons.



Backup



- $\gamma + \gamma$ separation criteria:
 - Two gammas (5GeV): varying distance
 - Efficiency definition: successful reconstruction of at least 2 neutral particles, both in 3.3GeV<E<6.6GeV
 - Removed events with γ -conversion before entering ECAL
 - Applied energy calibration
- $\gamma + \pi^+$ separation criteria:
 - 10GeV π + and 5GeV γ : varying distance
 - 3 T magnetic field
 - π + momentum measured by tracker
 - Efficiency definition: successful reconstruction of 3.3GeV<EN<6.6GeV, 9.9GeV<EC<10.1GeV
 - Removed events with γ/π + interactions before entering ECAL
 - Applied energy calibration

Sketch of ECAL in r-z plane







Performance studies: neutral pions with Arbor-PFA

- Reconstruction of π^0 in crystal ECAL: invariant mass and its resolution
 - Single π^0 's generated by the particle gun

Zhiyu Zhao (IHEP/SJTU)



Reconstruction: 3D

- Longitudinal linking for 3D cluster:
 - Cone-based clustering algorithm.
 - Get the very preliminary 3D structure, identify the cluster (MIP/EM/Hadron) with lateral

moment: $LAT = \frac{\sum_{i=3}^{N} E_i r_i^2}{\sum_{i=3}^{N} E_i r_i^2 + E_1 r_0^2 + E_2 r_0^2}$

• MIP: $\sqrt{LAT_X^2 + LAT_Y^2} < 0.05$

• EM:
$$0.05 < \sqrt{LAT_X^2 + LAT_Y^2} < 0.12$$

• Had:
$$\sqrt{LAT_X^2 + LAT_Y^2} > 0.12$$

PID Truth 10GeV	MIP	EM	Had
MIP	0.975	0.01	0
EM	0	0.99	0.01
Had	0.44	0.1	0.275





Measurements of the BGO energy resolution: setup

BGO Crystal:

- Lengths: 40/80/160mm
- Widths: 20/15/10mm
- Surfaces: polished/ground
- Tyvek / ESR wrapping





 4×4 mm² window for SiPM readoout

Photosensitive Device:

- SiPM & PMT
- SiPM: S13360-3050CS
 - 50 μ m pitch, 3 × 3mm², 3600 pixels
- PMT: R11065
 - 76mm (3"), gain: 5×10⁶





