LIDINE 2025

A Compact Optical Liquid Argon Facility at Roma Tre


Valerio D'Andrea¹, Giuseppe Salamanna^{1,2}, Hexi SHI (Kaku)^{1,2}, Diego Tagnani^{1,2}


- 1. Istituto Nazionale di Fisica Nucleare, Sez. RomaTre
- 2. Università degli Studi Roma Tre


Large Enriched Germanium Experiment

OLAF for LEGEND-1000

Germanium Experiment for Neutrinoless ββ Decay







- L-200 operating since 2022;
- up to 200 kg Ge;
- LAr as active shielding

- L-1000 baseline design ready;
- ~1000 kg Ge;
- LAr as active shielding



explore/qualify readout solutions for large assembly improve understanding of LAr as scintillation detector

LEGEND LAr talks:
Lucia Canonica
Luigi Pertoldi
Rosanna Deckert

OLAF Optical Liquid Argon Facility

setup at Roma3 featuring:

LAr: Ar 6.0 gas liquefied by LN₂

LAr scinti. photon by SiPM(VUV)

OLAF

Optical Liquid Argon Facility

benchmark/test:

LAr - SiPM - Front End - Trigger - DAQ readout chain for Legend-1000

setup at Roma3 featuring:


LAr: Ar 6.0 gas liquefied by LN₂

LAr scinti. photon by SiPM(VUV)


LAr optical properties that can be studied


- scintillation photon yield
- (effective) attenuation length
- time constants
- impurity effect (wo purification or doping)
- reflective surfaces

OLAF setup overview




OLAF setup overview





OLAF in preparation at INFN-Roma3


configuration for Ar liquefaction test




vacuum and cryogenic vessels, LN₂ dewar

LDetek_LD8000 Trace gas impurity analyser


LAr liquefaction with LN₂




LAr liquefaction with LN₂



Source and detector configuration - LLAMA case





Light source and trigger


<u>Am-241</u>

triple coincidence of 3 SiPMs near source


LED (green)


external pulser for power supply


OLAF: SiPM and light source mounting structures



OLAF: SiPM and light source mounting structures







OLAF Am-241 source housing and trigger


√ Am-241 gamma source at Roma3:


Commercially sealed capsule

Gamma: mainly 60 keV

Activity: 200 kBq

Activity tolerance: ± 30 %





OLAF Am-241 source housing and trigger


√ Am-241 gamma source at Roma3:


Commercially sealed capsule

Gamma: mainly 60 keV

Activity: 200 kBq

Activity tolerance: ± 30 %






LED source

- single photon calibration
- timing calibration
- basic functionality validation

SiPM TRG: in sync with LED power modulator




5 mm


Table 1a. Product performance of LUXEON Z Color Line at 500mA, T_i=25°C.

COLOR	DOMINANT [1] or PEAK WAVELENGTH [2] (nm)		LUMINOUS FLUX [1] (lm) or RADIOMETRIC POWER [3] (mW)		PART
	MINIMUM	MAXIMUM	MINIMUM	TYPICAL	NUMBER
Deep Red	650	670	250	350	LXZ1-PA01
Red	620	645	40	52	LXZ1-PD01
Red-Orange	610	620	56	65	LXZ1-PH01
Green	520	540	80	118	LXZ1-PM01
Cyan	490	510	48	82	LXZ1-PE01
Blue	460	480	24	38	LXZ1-PB01
Royal Blue	440	460	500	575	LXZ1-PR01

Trigger logic with V2740 digitizer Firmware



Trigger logic with V2740 digitizer Firmware



SiPM readout: Front-End, trigger, and DAQ

2023 JINST 18 P09007



LEGEND-200 FE boards x 2 & controller board

12 x 2 **RJ45** diff. outputs D. Tagnani

FE output adaptor

12 x 2 diff. outputs



CAEN V2740 digitizer



up to 24 SiPM channels

- connection to SiPM
- controller interface
 - -> validation to follow

FPGA FW for TRG logic MIDAS based DAQ


in preparation

DAQ implementation and first FW test

CAEN Dig V2740 FELib OEM FW for PSD




MIDAS package of TRIUMF Linux, C++ library; Online Data Base (ODB) for control and logging; GUI control from browser




TRG logic and run-time settings of digitizer —> set from GUI

64 CH 8 ns/S

DAQ basic validation: timing structure





for online PSD application



Ongoing prototype under test


assembly







Ongoing prototype under test

assembly







Geant4 simulation preparation

Ground work for G4 MC


"minimalist" setup implemented:


- LAr vessel, LAr volume
- tentative positions of 12 SiPMs (as Si boxes)
- geantino runs for geometry acceptance
- Optical Properties & realistic geometry under construction

number of hits @ SiPMs

1400 | Student |

geantino hits @ SiPM





Summary

Optical Liquid Argon Facility @ Roma Tre

- As an R&D test bench for Legend-1000 LAr detector
- Cryogenic setup constructed
- LAr scintillation direct readout with Hamamatsu VUV4 SiPMs

On-going activities:

- Argon liquefaction procedure establishment
- SiPM frame structure prototype production
- SiPM readout full chain preparation:
 SiPM -> Legend FE boards -> CAEN V2740 digitizer
- FPGA based trigger logics (using digitizer FW)
- DAQ: TRIUMF MIDAS based system

notes

Acknowledgement

On design thanks to the LLAMA team at TUM:

Patrick Krause, Laszlo Papp, Stefan Schönert, and Mario Schwarz

Christoph Vogl (working on SCARF, Xenon doping, and LAr purification)

Possible joint effort with Llama-SCARF


measurement with: reflective surface, doping, impurity ...

Mechanical workshop team at Roma3: Gianfranco Paruzza and Co.


HAMMER team at Roma1 for 3D printing

Special thanks to: Krzysztof Szczepaniec @ LNGS


External TRG readout test



Triple coincidence self-trigger readout test



V2740 DPP-PSD (charge) integration gate length



G4 geometry preparation

Current environment

MacOS 10.13.6 Geant4.10.2

