



NED 2026



LGAD用于BNCT中子束流测量

Study of LGAD sensor for BNCT Neutron Beam Measurement

Mengzhao Li (李梦朝), Bin Zhou, Tianjiao Liang, Mei Zhao, Zhijun Liang

Spallation Neutron Source Science Center
Institute of High Energy Physics, CAS

11, Aug. 2026

email: mzli@ihep.ac.cn

BNCT简介



硼中子俘获治疗(BNCT): 一种治疗癌症的新利器

- ① 硼药(^{10}B)在癌细胞富集
- ② 中子束流照射
- ③ 俘获反应产生 α 粒子和 ^7Li
- ④ 选择性杀死癌细胞

加速器BNCT装置:

- 日本8家, 至少2个已商业化并纳入医保
- 中国建成和在建有13家, 规划中约15家



东莞人民医院: 国产RFQ-BNCT一期临床完成

厦门弘爱医院: I期注册临床试验全部入组完成

青岛齐鲁医院: 回旋加速器BNCT启动临床试验

湖州市中心医院: 临床试验获批

.....



BNCT简介

BNCT治疗束的精确测量对疗效和患者安全至关重要

BNCT中子束参数:

- 超热中子 (0.5 eV–10 keV)
- 中子通量 $> 1 \times 10^9 \text{ cm}^{-2} \cdot \text{s}^{-1}$
- γ 本底注量率约 $10^7 \text{ cm}^{-2} \cdot \text{s}^{-1}$



束测要求:

- n/ γ 甄别
- 抗辐照能力
- 高计数率
- 空间分布扫描


BNCT中子束流探测器

- ^3He 管、金刚石
- 锂玻璃、LiCAF晶体


以进口探测器为主，需求量大、价格高、售后困难




低增益雪崩探测器(LGAD)简介



LGAD (Low-Gain Avalanche Diode)



LGAD electric field distribution



Single pixel IHEP-LGAD sensor


IHEP (Institute of High Energy Physics) LGAD sensor design :

- Active thickness 30-50 μm $\rightarrow$ Short carrier drift time $\rightarrow$ Fast signal rise time $\rightarrow$ Ultra-fast pulse ~ 3 ns
- Gain layer under PN junction $\rightarrow$ Low gain 10-50 $\rightarrow$ Lower shot noise $\rightarrow$ Higher SNR $\rightarrow$ Time resolution ~ 30 ps
- Shallow carbon doping process $\rightarrow$ Protect the gain layer $\rightarrow$ Improve radiation hardness $\sim 10^{16}$ $n_{\text{eq}}/\text{cm}^2$


LGAD is a promising domestic option for BNCT beam measurement.

低增益雪崩探测器(LGAD)简介


Timeline of IHEP-LGAD R&D




IHEP-NDL (2019)




IHEP-IMEv1 (2020.9)




IHEP-IMEv2 (2021.6)




IHEP-IMEv3 (2022.5)



Pre-production for ATLAS (2023.7)



Mass production (2024.7) ongoing



IHEP-LGAD has been deployed in the CERN ATLAS HGTD, and is under active exploration for BNCT.

LGAD+转换层探测热中子




LGADs are designed for charged particles, how to detect thermal neutrons ?


➤ Thermal neutron convertor materials

- ✓ Coating neutron convertor on LGAD surface, such as Lithium-6 (${}^6\text{LiF}$), Boron-10 (${}^{10}\text{B}$)
- ✓ The reaction between thermal neutrons and Li-6 produces Triton and α particles. Triton or α enter the LGAD sensitive region and generate neutron signals.
- ✓ Alpha and Triton will deposit all energy in 50 μm LGAD, and its signal is very large
- ✓ The electrons generated by gamma can only deposit a small amount of energy, and its signal is very small

LGAD coated with convertor ${}^6\text{LiF}$



Neutrons react with Li-6 and B-10




基于LGAD的BNCT中子探测器模块



LGAD based BNCT beam detector

Sensor parameter: IHEP-LGAD


- Sensitive area: 2.6 mm × 2.6 mm
- Sensitive thickness: 50 μm
- ^6LiF thickness: 2 μm



IHEP-LGAD sensor

Detector module parameter




- Size: 50 mm × 10 mm × 9 mm
- High voltage and signal share one SMA interface



LGAD detector module

基于LGAD的BNCT中子探测器模块

Alpha source test in the Lab Curium (Cm244) source: 5.8MeV α



- LGAD has a fast signal with a width of about 4 ns
- In the PHS of alpha, the pulse height of 5.8 MeV α is 176.2 mV, and the energy resolution is about 1.7% (sigma/mean)



LGAD模块在BNCT的测试

BNCT Equipment @ Dongguan



- | | |
|-----------------------------|---------------------------|
| ① Power source | ④ Treatment Room 1 |
| ② RFQ Accelerator | ⑤ Treatment Room 2 |
| ③ Power and control cabinet | ⑥ Vertical treatment room |

BNCT Equipment @ Dongguan People's Hospital

- This BNCT equipment was built by CASBNCT company
- Technology transfer from IHEP and CSNS
- Proton beam provided by 2.8 MeV RFQ linear accelerator
- The proton beam hits the lithium target to produce high-flux neutron beam

Neutron beam parameter

- Frequency: 200 Hz
- Duty ratio: 50 %




LGAD模块在BNCT的测试


Test setup at BNCT treatment room

LGAD detector module:


- LGAD with ${}^6\text{LiF}$ 2 μm




The neutron beam passes through the PMMA (polyethylene) moderator and hits the LGAD detector



BNCT treatment room
@ Dongguan People's Hospital




LGAD detector modules




LGAD模块在BNCT的测试

The properties of LGAD signals




LGAD signals in one beam cycle

- One beam cycle: 5 ms
- Neutron pulse width: 2.5 ms (duty ratio 50%)
- These signals contain neutrons and background




The LGAD signals

- Signal width ~ 4 ns
- Rise time ~ 1 ns
- Noise RMS < 2 mV
- SNR of neutron signal ~ 50




LGAD模块在BNCT的测试

PHSs of LGAD with ^6LiF



LGAD PHSs at beam ON & OFF

Neutron flux measurement



Neutron flux at different water depths

❑ Beam ON, the LGAD PHS consists of 4 parts:

- The preliminary analysis should be triton, alpha particle, BG. γ , and elastic scattering (ES) of n and Li

❑ Beam OFF, the LGAD PHS has only background γ


❑ LGAD measurements align with gold wire activation, validating LGAD for BNCT flux monitoring.



- ❑ LGAD is a new ultra fast silicon detector with strong radiation hardness
- ❑ The LGAD developed by IHEP has been used in the CERN ATLAS experiment
- ❑ We have developed a LGAD based neutron beam detector for BNCT
- ❑ Preliminary verification of the feasibility of LGAD for neutron measurement in BNCT

Next step

- Detection efficiency calibration
- Long time stability test (BNCT and CSNS beam)
- More BNCT beam test (fluence, profile, ect)



LGAD detector module

欢迎合作！ 共同拓展LGAD的应用研究。

THANK YOU!

高薪!

中国散裂中子源大气中子辐照谱仪公开招聘两名博士后:

岗位	岗位职责	岗位要求
岗位1: 高能中子能谱测量及方法研究	聚焦 高能中子探测 ，发展高能区中子测量新方法；开展新型中子探测器研发及多场景大气中子能谱测量；并参加谱仪日常工作	取得博士学位3年以内；核物理、粒子物理、辐射探测相关专业；具备辐射探测器研发、中子能谱解谱、蒙特卡罗输运仿真经验
岗位2: 先进半导体器件与辐照效应研究	面向 先进工艺半导体器件及设备 ，开展大气中子辐照效应实验及机理研究；并参加谱仪日常工作。	取得博士学位3年以内；微电子、半导体物理、辐射效应相关专业；熟练使用TCAD及蒙卡程序，具有单粒子效应实验经验