



中国科学院大学
University of Chinese Academy of Sciences

Study of Thinner Thick Gaseous Electron Multiplier

Zhang Yuning, Q. Liu, H. B. Liu, Y. H. Zheng,
Y. G. Xie, W. Q. Huang, B. L. Wang, S. Chen,
W. T. Luo

5th Micro-Pattern Gas Detector Workshop



Outline

- ◆ Introduction
- ◆ Spatial resolution
- ◆ Sealed chamber
- ◆ Summary

Thick Gaseous Electron Multiplier (THGEM)



- ◆ THGEMs
 - Thickness $t = 0.4 \sim 3 \text{ mm}$
 - Hole diameter $d = 0.3 \sim 1 \text{ mm}$
 - Pitch $a = 0.7 \sim 7 \text{ mm}$
- ◆ Robust
- ◆ Can be cascaded for higher gain
- ◆ Effective single-photon detection in cascade + photocathode
- ◆ Few-ns RMS time resolution
- ◆ Cryogenic operation: OK
- ◆ Sub-mm spatial resolution (for $t=0.4 \text{ mm}$, $d=0.5 \text{ mm}$, $a=1 \text{ mm}$, spatial resolution is 0.7 mm)

*L. Arazi, DARWIN meeting WIS Jan 2015




Thinner-THGEM

- ◆ Thinner-THGEM: $t = 0.2$ mm, $d = 0.2$ mm, $a = 0.5$ mm, rim = $5 \sim 20 \mu\text{m}$.
- ◆ Thinner-THGEM chamber: active area of $5 \times 5 \text{ cm}^2$.



- ◆ Advantages:
 - Under the conditions of obtaining the same gain, the operating voltage is **lower**.
 - It is easy to **curve**. So it can be used for one-dimensional X-ray diffraction imaging.
 - **Better spatial resolution.**

Spatial resolution





$$\sigma^2 = \frac{\sum (x_i - \bar{x})^2}{N}$$



Impact factors

1. Structure of cathode
2. E_{drift}
3. THGEM's parameters
4. E_{Ind} ,
5. Anode straps' size
6. Electronics

Experiment setup





Test procedure





- ◆ THGEM:
150 μ m/300 μ m
- ◆ Anode:
150 μ m/300 μ m
- ◆ We deem x_0 is the middle of the slit (the width of slit is 20 μ m)

Centers method



- ◆ Event by event.
- ◆ Digit signal
- ◆ Gastone64 read out.
- ◆ $\sigma_{\min} = 0.5712 \times 300 \mu\text{m} = 171.36 \mu\text{m}$
- ◆ $\sigma_{\max} = 0.6917 \times 300 \mu\text{m} = 207.51 \mu\text{m}$

Photo Multiplier Tube (PMT)



◆ PMT



- High Gain
- Excellent time resolution (ps)
- Output channel limited
- Magnetic field deflected
- Expensive

◆ GPM

- large areas, flat geometry
- operation in magnetic fields
- sensitivity to single photons
- spectral range from UV to visible
- fast (ns range)
- high localization accuracy (sub-mm range)



Seal

- ◆ Visible light
- ◆ Semi-transparent photocathodes
- ◆ K-Cs-Sb photocathodes are very chemically reactive and decay promptly. Therefore, detector comprising bialkali photocathodes must operate in a sealed chamber.
- ◆ Astrophysics (Atmospheric Cherenkov), medical applications, ...




Experiment setup

- ◆ We designed a small sealed chamber to study various types of anode, outgas of THGEM foils and so on.
- ◆ The front end uses a quartz of $\varnothing 80 \times 2 \text{ mm}$ as the transparent window. Use the aciculiform valves to seal.
- ◆ Chase leaks with He. ($10^{-10} \text{ Pa} \times \text{m}^3/\text{s}$)
- ◆ Ar (97%) + iC_4H_{10} (3%)





X-ray Spectrum

histogram of MCA




histogram of MCA




- ◆ Cu spectrum by gas flow THGEM detector: full-energy peak and escape peak.
- ◆ Cu spectrum by CdTe detector: 8.04 keV (K_{α} , 80%) and 8.9 keV (K_{β} , 20%).

X-ray Spectrum


histogram of MCA



histogram of MCA





SiO2(2mm)



- ◆ High energy Bremsstrahlung photons which come from the X-ray tube, can be transmitted through the quartz window, and hit the copper layer of the THGEM foil.

Gain



- ◆ The bigger the pressure of the chamber, the gain is smaller.
- ◆ Long-term stability. The data were obtained at a pressure of 60 kPa.
- ◆ Pressure remain unchanged.
- ◆ The gain decreases ~30% in 20 days.

Current state



Summary and next plan

- ◆ Spatial resolution: $\sigma = 171.36 \mu\text{m}$ (Centers method).
- ◆ Preliminary study sealed chamber.
- ◆ We will further study center of gravity method.
- ◆ We will further study the technology of coating alkali-antimonides.

Thank you!