



Ultralow background LS and related technologies

Grzegorz Zuzel

Jagiellonian University, Cracow, Poland

Outline



- Ultra-pure LS: BOREXINO
- Internal and external background: mitigation techniques
- LS purification
- Conclusions



BOREXINO at LNGS




BOREXINO

Bcg mitigation

LS purification

Conclusions





BOREXINO design




BOREXINO

Bcg mitigation

LS purification

Conclusions





BOREXINO radio-purity

In a nutshell: the cleanest detector ever built



BOREXINO

Bcg mitigation

LS purification



| Isotope Specification for LS | | Achieved after filling (2007 - 2010) |
|----------------------------------|-----------------------------|--|
| 238 U | $\leq 10^{-16} \text{ g/g}$ | $(5.3 \pm 0.5) \cdot 10^{-18} \text{ g/g}$ |
| ²³² Th | $\leq 10^{-16} \text{ g/g}$ | $(3.8 \pm 0.8) \cdot 10^{-18} \text{ g/g}$ |
| ¹⁴ C/ ¹² C | ≤ 10 ⁻¹⁸ | $(2.69 \pm 0.06) \cdot 10^{-18} \text{ g/g}$ |
| $^{40}\mathrm{K}$ | $\leq 10^{-18} \text{ g/g}$ | $\leq 0.4 \cdot 10^{-18} \text{ g/g}$ |
| ⁸⁵ Kr | ≤ 1 cpd/100 t | $(30 \pm 5) \text{ cpd/}100 \text{ t}$ |
| ³⁹ Ar | ≤ 1 cpd/100 t | << ⁸⁵ Kr |
| ²¹⁰ Po | not specified | ~ (70) 1 dpd/100 t |
| ²¹⁰ Bi | not specified | (20) 70 dpd/100 t |

$$A_{Bx} = \sum (bcg \ ev.) \sim 30 \frac{cpd}{100 \ t} \sim \mathbf{10^{-9}} \frac{Bq}{kg}$$

$$A_{water} \sim \mathbf{10} \frac{Bq}{kg}$$

$$\rightarrow activity \ reduction \ factor$$

$$f \sim \mathbf{10^{10}}$$





BOREXINO

Bcg mitigation

LS purification



- Graded shielding: traveling inward to the center, each component is protected from external radiation by the preceding one
- The radio-purity level is increasing towards the center
- Active (definition of FV, Čerenkov veto) and passive (PC buffer, water) suppression of external radiation
- Careful selection of construction materials and detector components with respect to content of radioactive isotopes, ²²²Rn emanation and permeability
- Preventing surface contamination
- Application of appropriate purification (liquids, gases) and cleaning techniques

BOREXINO design




BOREXINO

Bcg mitigation

LS purification

Conclusions





CTF – testing the scintillator



BOREXINO

Bcg mitigation

LS purification

Conclusions

Nn[,]16



CFT – 1: 1995

CTF - 2:2000

CTF - 3: 2001-2003

- Phys. Lett. B, 422, 349 (1998)
- Astrop. Phys. 8(3), 141 (1998)
- NIM A406, 411 (1998)
- Physics Letters B 525, 29 (2002)
- Physics Letters B 563, 23 (2003)
- Physics Letters B 563, 37 (2003)
- JETP Lett. 78 No 5, 261 (2003)
- Eur. Phys. J. C 37, 421 (2004)
- Eur. Phys. J. C 47, 21 (2006)
- Phys. Rev. C 74, 045805 (2006)

$$\frac{^{14}C}{^{12}C} \sim 10^{-18}$$
 $C_{U/Th} \sim 4 \times 10^{-16} \ g/g$

HPGe spectroscopy



BOREXINO

Bcg mitigation

LS purification

Conclusions



GeMPIs at GS (3800 m w.e.)


- GeMPI I operational since 1997 (MPIK)
- GeMPI II built in 2004 (MCavern)
- GeMPI III/IV constructed in 2007-2012 (MPIK/LNGS)
- World's most sensitive spectrometers

GeMPI I:

- Crystal: 2.2 kg, $\varepsilon_{\rm r} = 102 \%$
- Bcg. Index (0.1-2.7 MeV): 6840 cts/kg/year
- Sample chamber: 15 l

Sensitivity for U/Th: ~10 μBq/kg

Appl. Rad. Isot., 53 (2000) 191 Astrop. Phys. 18 (2002) 1



Survey different materials

| X | Septem 1 |
|---|----------|
| | X |

BOREXINO

Bcg mitigation

LS purification

Conclusions



| | Sample | Description | ²²⁶ Ra [mBq/kg] | ²²⁸ Th [mBq/kg] | ⁴⁰ K [mBq/kg] |
|--|--------------------|------------------|-------------------------------|-------------------------------|----------------------------|
| | a | AISI304L: SSS | 4.6 ± 0.9 | 11.4 ± 1.1 | < 14 |
| | Stainless Steel | SS for pipes | < 14 | < 10 | < 34 |
| | Steel | SS for flanges | 6.2 ± 1.2 | 6.5 ± 1.6 | < 13 |
| | PMTs | Dynodes | < 280 | 450 ± 163 | < 240 |
| | inner | Ceramic plates | 170 ± 50 | 310 ± 60 | 960 ± 450 |
| | parts | Al for dynodes | 1190 ± 100 | 980 ± 80 | 2800 ± 600 |
| | PMTs | Mu metal | 57 ± 20 | < 27 | < 180 |
| | ancillary parts | Volt. div. board | 170 ± 60 | 80 ± 40 | 770 ± 360 |
| | | Voltage divider | 680 ± 30 | 320 ± 20 | 3200 ± 320 |
| | Glass | Sand for glass | 40 ± 3 | < 3.1 | < 25 |
| | | ETL LB glass | 820 ± 230 | 130 ± 12 | 500 ± 120 |
| | | Base glass | 520 ± 90 | 410 ± 90 | $(2.2 \pm 0.6) \cdot 10^5$ |
| | | | | | |

Astroparticle Physics 18 (2002) 1

²²²Rn emanation




BOREXINO

Bcg mitigation

LS purification

Conclusions

Nn[,]16



²²²Rn emanation



BOREXINO

Bcg mitigation

LS purification

Conclusions





Blanks:

 $201 \rightarrow 50 \mu Bq$

 $801 \rightarrow 80 \mu Bq$

Absolute sensitivity ~100 μBq [50 atoms]

Appl. Rad. Isot. 53 (2000) 371

²²²Rn emanation: examples



BOREXINO

Bcg mitigation

LS purification

Conclusions



| System | Sample | Description | ²²² Rn em. rate |
|--|---|-------------------------------------|-------------------------------|
| | SS vessel TK1 | 114 m ³ | < 60 mBq |
| PC storage area | SS vessel TK2 | 114 m^3 | $(45 \pm 8) \text{ mBq}$ |
| | SS vessel TK3 | 114 m^3 | $(24 \pm 5) \text{ mBq}$ |
| N_2 | Electrical heater | | $(0.92 \pm 0.29) \text{ mBq}$ |
| distribution line LS purification plant | Particle Filter | | $(0.34 \pm 0.13) \text{ mBq}$ |
| | 1.5" distrib. line | ~ 100 m long | $(0.47 \pm 0.13) \text{ mBq}$ |
| | SS package | 25 m^2 | < 0.12 mBq |
| | H ₂ O extraction column + 24 SS packages | $0.6 \text{ m}^3 / 608 \text{ m}^2$ | $(4.83 \pm 0.70) \text{ mBq}$ |
| | N ₂ sparging column + 26 SS packages | $0.2 \text{ m}^3 / 280 \text{ m}^2$ | $(1.78 \pm 0.21) \text{ mBq}$ |

Over 1000 entries in the DB!

Astroparticle Physics 18 (2002) 1 LRT 2004 proceedings, p. 141 – 149 Int. J. Mod. Phys. A29 (2014) 1442009

²²²Rn/²²⁶Ra in water



BOREXINO

Bcg mitigation

LS purification

Conclusions





International

STRAW: System for the ²²²Rn and ²²⁶Ra Assay of Water

- Placed at the BOREXINO water plant
- 222Rn extraction from 350 liters
- 222Rn and 226Ra measurements possible

²²²Rn detection limit: ~0.1 mBq/m³ ²²⁶Ra detection limit: ~0.8 mBq/m³

Nucl. Instr. Meth. A 497 (2003) 407

5 Beijing, China

| H ₂ O flow [m ³ /h] | HPN ₂ flow [kg/h] | C_{Rn} [mBq/m ³] | ${ m C_{Ra}} \ [{ m mBq/m^3}]$ |
|---|---------------------------------|--------------------------------|--------------------------------|
| 2 30 | | 704 ± 7 | 1.2 ± 0.5 |
| 1 20 | | 247 ± 6 | 3.8 ± 0.7 |
| 1 30 | | 186 ± 5 | 2.0 ± 0.6 |
| Loop mode | | 3.0 ± 0.4 | 1.3 ± 0.9 |

BOREXINO design




BOREXINO

Bcg mitigation

LS purification

Conclusions





²²²Rn diffusion




BOREXINO

Bcg mitigation

LS purification

Conclusions





Sensitivity: D ~ 10^{-13} cm²/s d_e ~ 2 μ m

²²²Rn diffusion



Results obtained for the 0.018 mm thick C38F film (BOREXINO)

| RH standard salt | RH in gas phase (%) | Water amount in nylon, M (%) | Diffusion coefficient, $D \text{ (cm}^2/\text{s)}$ | Solubility, S |
|--|------------------------|--------------------------------|--|------------------|
| Mg(ClO ₄) ₂ | ~0 | ~0 | $(2.1\pm0.4)\times10^{-12}$ | 4.5 ± 0.7 |
| $H_3PO_4 \cdot \frac{1}{2}H_2O$ | 9 ± 1 | 0.72 ± 0.04 | $(2.3\pm0.3)\times10^{-12}$ | 2.5 ± 0.3 |
| $\text{LiCl}_2 \cdot \overline{\text{H}_2\text{O}}$ | 12 ± 1 | 0.87 ± 0.04 | $(2.2\pm0.3)\times10^{-12}$ | 2.2 ± 0.3 |
| CaCl ₂ · 6H ₂ O | 32 ± 2 | 2.09 ± 0.04 | $(4.3\pm0.5)\times10^{-12}$ | 1.8 ± 0.2 |
| $Na_2Cr_2O_7 \cdot 2H_2O$ | 52 ± 2 | 3.74 ± 0.05 | $(1.9\pm0.3)\times10^{-11}$ | 1.4 ± 0.2 |
| $Na_2S_2O_3 \cdot 5H_2O$ | 76 ± 2 | 6.35 ± 0.05 | $(6.5\pm0.9)\times10^{-11}$ | 1.5 ± 0.2 |
| K ₂ CrO ₄ | 88 ± 3 | 7.60 ± 0.05 | $(1.3\pm0.2)\times10^{-10}$ | 1.5 ± 0.2 |
| Na ₂ SO ₄ · 10H ₂ O | 93 ± 3 | 9.12 ± 0.07 | $(3.3\pm0.4)\times10^{-10}$ | 1.0 ± 0.1 |
| H ₂ O vapors | 100 ± 3 | 10.14 ± 0.09 | $(1.3\pm0.2)\times10^{-9}$ | 0.7 ± 0.1 |

BOREXINO

Bcg mitigation

LS purification

Conclusions



There is 3 orders of magnitude difference between the diffusion in the dry and in the foil saturated with water!

Nucl. Instr. Meth. A 449 (2000) 158

Nucl. Instr. Meth. A 524 (2004) 355

$$d_e = \sqrt{\frac{D}{\lambda}}$$
 $d_e^d = 7 \mu m$ $d_e^w = 270 \mu m$

²²⁶Ra in/on BOREXINO nylon



BOREXINO

Bcg mitigation

LS purification

Conclusions





1 ppt U required (~12 μBq/kg for ²²⁶Ra)

 $D_{dry} = 2x10^{-12} \text{ cm}^2/\text{s} \ (d_{dry} = 7 \text{ } \mu\text{m})$ $D_{wet} = 1x10^{-9} \text{ cm}^2/\text{s} \ (d_{wet} = 270 \text{ } \mu\text{m})$

 $A_{dry} = A_{sf} + 0.14 \cdot A_{bulk}$ $A_{wet} = A_{sf} + A_{bulk}$

Separation of the bulk and surface ²²⁶Ra conc. was possible through ²²²Rn emanation

Very sensitive technique: $(C_{Ra} \sim 10 \mu Bq/kg)$

Bx IV foil: bulk $\leq 15 \mu Bq/kg$ surface $\leq 0.8 \mu Bq/m^2$ total = $(16 \pm 4) \mu Bq/kg (1.2 ppt U eqiv.)$

NIM A 498 (2003) 240

Construction of nylon vessels



BOREXINO

Bcg mitigation

LS purification

Conclusions





Princeton clean room class 100 with 222 Rn-reduced air (VSA filter): $C_{Rn} \sim 1 \text{ Bq/m}^3$

A. Pocar, PhD Thesis (2003)

Inflation of vessels in SSS




BOREXINO

Bcg mitigation

LS purification

Conclusions





The nylon vessels were inflated in the sphere with synthetic air: $C_{Rn} < 100 \mu Bq/m^3$

Int. J. Mod. Phys. A29 (2014) 1442009





BOREXINO

Bcg mitigation

LS purification

Conclusions



Effective scintillator purification processes


- Efficient removal of radioactive elements
 - distillation
 - water extraction
 - nitrogen stripping
- Development of cleaning methods to remove radioactivity from equipment surfaces
- Purification of scintillator before and after filling detector



- **BOREXINO**
- Bcg mitigation
- LS purification
- Conclusions
- Nn[,]16

- Vessels inflated with synthetic air
- Vessels filled with high-purity deionized water
- Scintillator distilled and stripped with LAK nitrogen during filling
- Detector filled with scintillator from top, while draining water from bottom

J. Benziger et al. / Nuclear Instruments and Methods in Physics Research A 587 (2008) 277-291



BOREXINO design




BOREXINO

Bcg mitigation

LS purification

Conclusions





²²²Rn in gases: MoREx




BOREXINO

Bcg mitigation

LS purification





Ar and Kr in nitrogen




BOREXINO

Bcg mitigation

LS purification





BOREXINO nitrogen



BOREXINO

Bcg mitigation

LS purification

Conclusions



Regular Purity Nitrogen:

- Technical 4.0 quality, not purified
- Production rate up to 100 m3/h (STP)
- 222 Rn (30 70) µBq/m³ Ar ~ 10ppm, Kr ~ 30 ppt

High Purity Nitrogen:

- ²²²Rn adsorption on charcoal (LTA)
- Achieved concentration (0.30 ± 0.09) μBq/m³
- Production rate up to 100 m³/h (STP)
- Ar and Kr not removed



LAK (Low Ar and Kr) Nitrogen:

- Spec. Ar < 0.4 ppm, Kr < 0.2 ppt 222Rn < 7 μBg/m³
- Purification by adsorption on different materials extensively studied (successfully!)
- Cooperation with companies on the nitrogen survey
- Tests of the nitrogen delivery chain

Nitrogen survey

| Nitrogen sample | C _{Ar} [ppm] | C _{Kr} [ppt] | |
|--------------------|-----------------------|-----------------------|--|
| MESSER (4.0) | 200 ± 30 | 1680 ± 240 | |
| Air Liquide (4.0) | 11.0 ± 1.3 | 40 ± 5 | |
| Linde AG, (7.0) | 0.031 ± 0.004 | 2.9 ± 0.4 | |
| SOL (6.0) | 0.0063 ± 0.0006 | 0.04 ± 0.01 | |
| Westfalen AG (6.0) | 0.00050 ± 0.00008 | 0.06 ± 0.02 | |
| Goal (BOREXINO) | < 0.4 | < 0.2 | |

Tests of the delivery chains



| Supplier/setup | $C_{Rn} \left[\mu B q/m^3 \right]$ | C _{Ar} [ppm] | C _{Kr} [ppt] |
|---|-------------------------------------|-----------------------|-----------------------|
| Linde AG, 3-m ³ movable tank | 1.2 | 0.018 | 0.06 |
| SOL, 16-m ³ tank | 8 | 0.012 | 0.02 |

BOREXINO LAK nitrogen



BOREXINO

Bcg mitigation

LS purification

Conclusions





LS re-purification



BOREXINO

Bcg mitigation

LS purification



- From 2007 to 2010 the ²¹⁰Bi background increased
 - Scintillator refills to compensate the vessel leak / others ?
 - ²¹⁰Bi rate: 20 → 70 cpd/100 t
 - 85Kr rate: 30 cpd/100 t (constant)
- To reduce background, scintillator was re-purified using two processes:
 - Water extraction to remove ²¹⁰Pb (²¹⁰Bi)
 - Nitrogen stripping to remove ⁸⁵Kr and other volatiles
- Six purification cycles were performed, each took ~1 month
 - Each operation processed all the scintillator in the detector once (320 m³).
 - Data were acquired to evaluate backgrounds after each operation.
 - Not all six operations were successful.

BOREXINO radio-purity



BOREXINO

Bcg mitigation

LS purification



| Isotope | Specification for LS | Achieved after filling (2007 - 2010) | After aditional purification |
|----------------------------------|------------------------------------|--|--|
| ^{238}U | $\leq 10^{-16} \text{ g/g}$ | $(5.3 \pm 0.5) \cdot 10^{-18} \text{ g/g}$ | $< 0.8 \cdot 10^{-19} \text{ g/g}$ |
| ²³² Th | $\leq 10^{-16} \text{ g/g}$ | $(3.8 \pm 0.8) \cdot 10^{-18} \text{ g/g}$ | $< 1.2 \cdot 10^{-18} \text{ g/g}$ |
| ¹⁴ C/ ¹² C | ≤ 10 ⁻¹⁸ | $(2.69 \pm 0.06) \cdot 10^{-18} \text{ g/g}$ | unchanged |
| ⁴⁰ K | $\leq 10^{-18} \text{ g/g}$ | $\leq 0.4 \cdot 10^{-18} \text{ g/g}$ | unchanged |
| ⁸⁵ Kr | $\leq 1 \text{ cpd/}100 \text{ t}$ | $(30 \pm 5) \text{ cpd/}100 \text{ t}$ | ≤ 5 cpd/100 t |
| ³⁹ Ar | ≤ 1 cpd/100 t | << ⁸⁵ Kr | << 85Kr |
| ²¹⁰ Po | not specified | ~ (70) 1 dpd/100 t | unchanged |
| ²¹⁰ Bi | not specified | (20) 70 dpd/100 t | $(20 \pm 5) \text{ cpd/}100 \text{ t}$ |

Topics not discussed



BOREXINO

Bcg mitigation

LS purification



- ²²⁶Ra adsorption on the nylon vessel (J. Rad. Nuc. Chem. 296 (2013) 639)
- ²²²Rn-daughters deposition on the nylon vessel (E. Harding, Princeton)
- Tests of ²¹⁰Pb removal from PC (J. Rad. Nuc. Chem. 296 (2013) 639)
- Online ²²²Rn monitoring with an electrostatic detector (NIM A 460 (2001) 272)
- Adsorption of noble gases on various porous materials (B. Freudiger, PhD Thesis (2003))
- NAA and ICP-MS measurements (Astroparticle Physics 18 (2002) 1)

Conclusions



BOREXINO

Bcg mitigation

LS purification

Conclusions



- BOREXINO has achieved an unprecedented background level in the liquid scintillator
- Strict quality control program including the assay of all components of the detector during its construction
- +10 years of R&D, many people/institutions involved
- Several detectors and experimental methods were developed allowing measurements even at a single atom level.
- Most of the developed techniques are world-wide most sensitive (Ge spectroscopy, ²²²Rn detection, ²²²Rn diffusion) and are applied in next-generation experiments (GERDA, XENON, DARKSIDE,...)