TIPP17 International Conference on Technology and Instrumentation in Particle Physics MAY 22-26, 2017 | BEIJING, PEOPLE'S REPUBLIC OF CHINA

Progress in Room-Temperature and Cryogenic Resistive THGEM-based detectors

Shikma Bressler Weizmann Institute of Science





THGEM Thick Gaseous Electron Multiplier

Operation principle

A. Breskin et al. NIM A598 (2009) 107

- Mechanically drilled FR4 plate Cu-clad on both sides
- Placed a few mm above the readout pads
- Voltage supplied to the electrodes
 - Creates a strong dipole field between the THGEM electrode amplification region
 - Drift field focuses the primary ionization into the holes
 - Induction field pull the avalanche induced electrons to the anode



Typical geometry ×10 GEM

 $a \sim 1~mm,\, d \sim 0.5~mm,\, h \sim 0.1~mm$



Thick Gaseous Electron Multiplier

Wish list

Random order

- Simple
- Robust
- Cost-effective
 - Production
 - Operation (etc. gas mixtures)
- Large-area
- Efficient
- Resolution: spatial, time, energy
- Broad dynamic range
- Discharge free
- Rate capabilities
- Industrially produced

Applications

Random order

- RICH devices M. Alexeev et al. 2012 JINST 7 C02014
- Cryogenic detectors for TPC in neutrino physics and rare-event searches

M. Resnati et al. 2011 J. Phys.: Conf. Ser. 308 012016 A. Bondar et al. 2011 JINST 6 P07008

- GPM for dark matter searches L. Arazi et al. Expected online publication in JINST: November 2015
- Medical imaging S. Duval et al. 2011 JINST 6 P04007
- Neutron/Gamma imaging in cargo inspection systems A. Breskin et al. 2012 JINST 7 C06008

I. Israelashvili et al. 2015 JINST 10 P03030

- Thin sampling elements for DHCAL S. Bressler et al. 2013 JINST 8 P07017
- And more...



Thick Gaseous Electron Multiplier

Wish list

Random order

- Simple
- Robust
- Cost-effective
 - Production
 - Operation (etc. gas mixtures)
- Large-area
- Efficient
- Resolution: spatial, time, energy
- Broad dynamic range
- Discharge free
- Rate capabilities
- Industrially produced

Hard - common for all MPGDs

• Efficient

- MIP detection gain of 10⁴
- Single photon detection gain of 10⁵
- Background event depositing 10²-10³ primary electron ⇒ reach the Ruether limit ⇒ Discharge



Thick Gaseous Electron Multiplier

Wish list

Random order

- Simple
- Robust
- Cost-effective
 - Production
 - Operation (etc. gas mixtures)
- Large-area
- Efficient
- Resolution: spatial, time, energy
- Broad dynamic range
- Discharge free
- Rate capabilities
- Industrially produced

Hard - common for all MPGDs

- Efficient
 - MIP detection gain of 10⁴
 - Single photon detection gain of 10⁵
- Background event depositing 10²-10³ primary electron ⇒ reach the Ruether limit ⇒ Discharge

Solution I - cascade

- Split multiplication between several stages
 - The filed in each stage is lower
 - The charge is divided between several holes
- More expensive
 - Not suitable for applications like DHCAL requiring thin elements



Thick Gaseous Electron Multiplier

Wish list

Random order

- Simple
- Robust
- Cost-effective
 - Production
 - Operation (etc. gas mixtures)
- Large-area
- Efficient
- Resolution: spatial, time, energy
- Broad dynamic range
- Discharge free
- Rate capabilities
- Industrially produced

Hard - common for all MPGDs

- Efficient
 - MIP detection gain of 10⁴
 - Single photon detection gain of 10⁵
- Background event depositing 10²-10³ primary electron ⇒ reach the Ruether limit ⇒ Discharge

Solution II - Resistive configuration

• Single stage

- Operated with cost-effective gas mixtures (Ar-based rather than Nebased)
- Targeting discharge self-quenching mechanism
- Rate capabilities 5% efficiency drop over 3 orders of incoming particle rates

THGEM Resistive configurations - background

THGEM

drift electrode

MIP

 $\mathsf{E}_{\mathsf{drift}}$

THGEM

- Cu-clad on both sides
- Operation with induction gap

THick-WELL (THWELL)

- Cu-clad on one side
- Electrode coupled directly to the anode

Resistive-WELL (RWELL)

- WELL coupled to the readout through resistive layer (10-20 MΩ)
- Readout separated from the readout by a thin insulating sheet
- Signal induced on the readout



JINST 8 (2013) C12012

7

Discharge energy



International Conference on Technology and Instrumentation in Particle Physics MAY 22-26, 2017 | BEIJING, PEOPLE'S REPUBLIC OF CHINA

Resistive configurations

- Discharge energy is quenched
- Discharge probability does not change



International Conference on Technology and Instrumentation in Particle Physics MAY 22-26, 2017 | BEIJING, PEOPLE'S REPUBLIC OF CHINA JINST 9 2014 P03005

S. Bressler



Resistive Plate WELL

• Single sided THick Gaseous Electron Multiplier (THGEM)

RPWELL

• Coupled to the readout anode through material of high bulk resistivity



Voltage supplied to the RPWELL electrode creates a strong field inside the holes - amplification region.

The signal is induced on the readout anode.

Discharge quenching is achieved by delaying the charges on the RP - local reduction of the field

9



JINST 8 (2013) P11004

RPWELL

Resistive Plate WELL

- Single sided THick Gaseous Electron Multiplier (THGEM)
- Coupled to the readout anode through material of high bulk resistivity

10

Main characteristics - room temperature

- High single-electron gain (>10⁵) with single element
- High detection efficiency (MIP ~99%)
- Moderate gain loss under high rate of incoming particle fluxes
- Discharge free operation
 - With Ne- and Ar-based gas mixtures
 - At muon and high rate pion beams
 - Under broad dynamic range of primary ionization



JINST 11 (2016) no. 01 P01005

JINST 11 (2016) no. 09 P09013

NIM A845 (2017) 2620265

International Conference on Technology and Instrumentation in Particle Physics MAY 22-26, 2017 | BEIJING, PEOPLE'S REPUBLIC OF CHINA

S. Bressler



Moving to $50 \times 50 \text{ cm}^2$

 $3 \times 3 \text{ cm}^2 \Rightarrow 10 \times 10 \text{ cm}^2 \Rightarrow 30 \times 30 \text{ cm}^2 \Rightarrow 50 \times 50 \text{ cm}^2$

Motivation

 Applications requiring cost-effective large area detectors
 with moderate spatial resolution

with moderate spatial resolution

- (Semi) Digital Hadronic Calorimeter (S)DHCAL
- Muon tomography for homeland security
- Volcanology and more



MAY 22-26, 2017 | BEIJING, PEOPLE'S REPUBLIC OF CHINA



RPWELL Results in the context of a (s)DHCAL

- 10×10 and 30×30 cm² prototypes
- 1×1 cm² pad readout
- 0.4 mm thick Semitron resistive plate
- 5 mm drift gap
- Ne/CH₄(5%), Ar/CO₂(7%), Ar/CH₄(5%)
- Muon and high rate pion beams
- APV25/SRS readout



Discharge free operation also under high pion flux





DD

Moving to $50 \times 50 \text{ cm}^2$

$3 \times 3 \text{ cm}^2 \Rightarrow 10 \times 10 \text{ cm}^2 \Rightarrow 30 \times 30 \text{ cm}^2 \Rightarrow 50 \times 50 \text{ cm}^2$



S. Bressler

WEIZMANN

Moving to $50 \times 50 \text{ cm}^2$

Assembly

- Strip readout (1 & 1.5 mm pitch)
- Silicate glass Resistive Plate (10⁹ Ωcm)
 [J.-b. Wang et al., NIM A621 151]
- Anode to glass coupling through graphite-epoxy layer (~MΩ)
- Gluing under vacuum







International Conference on Technology and Instrumentation in Particle Physics MAY 22-26, 2017 | BEIJING, PEOPLE'S REPUBLIC OF CHINA



Moving to $50 \times 50 \text{ cm}^2$

 $3 \times 3 \text{ cm}^2 \Rightarrow 10 \times 10 \text{ cm}^2 \Rightarrow 30 \times 30 \text{ cm}^2 \Rightarrow 50 \times 50 \text{ cm}^2$

Performance

- Studies conducted with X-Y scanner
 - Developed for the ATLAS sTGC
- Ag-target x-ray tube \Rightarrow 22 KeV photons
 - Penetrating the 3 mm FR4 cover (cathode)
- High intensity ($\sim 400 \text{ KH}_z/\text{mm}^2$)
- HV current monitoring synchronized with the source position





- Ar/CO₂(7%) ΔV_{RPWELL} = 1800 V Flux ~ 400 kH_z/mm² 30 sec/point
- Stable operation
- Some non-uniformities
 - under study



Moving to $50 \times 50 \text{ cm}^2$

Next steps

- In the lab
 - Investigate & solve the non-uniformities (correlated with the glass tiles)
- In-beam evaluation (soon)
 - Detection efficiency
 - Uniformity
 - Discharge probability
 - Position resolution
 - Energy resolution

- (S)DHCAL
- Project shared with LAPP and Democritus
- Build prototype
 - With Pad readout
 - With the MICROROC chip
 - Design and assembly procedure identical to the one used
- Test in the lab and in the beam
- •
- •
- Incorporate within a stack of ~20 sampling elements; MICROMEGAS and RPWELL



UV detection

Motivation

- Applications requiring large coverage of UV detection
- In room temperature and cryogenic temperatures
- Possibly with low ion back flow

Example

- RICH detectors (COMPASS 1 upgrade)
- Phoenix Hadron Blind Detector
- Cryogenic Gaseous Photon Multipliers (GPM) for
 - Multi-ton dark matter experiment
 - Multi-ton neutrino detectors
 - Homeland security, etc



UV detection

COMPASS RICH-I detector

- Hybrid detector: two-stage THGEM+MICROMEGAS
- Operation gain $\sim 10^5$ Exponential single-photon distribution
- Good time resolution (7 ns RMS)
- Experiment ongoing







DOI: 10.1109/NSSMIC.2015.7581949



Gaseous Photon Multiplier



Light&Charge amplifications





4π all Liquid TPC



Gaseous Photon Multiplier

Cascade structure

- To reduce discharge rate
- Possibility to reduce Ion Back Flow

Imaging capabilities

- Pixelated readout
- Low noise







Imaging at Room Temperature



Single event image (H2 lamp); Center of gravity (COG) calculated for each event

2D histogram of COGs of events comprising 160 PEs

GPM

30

50



GPM Imaging at Cryogenic Temperature



S. Bressler



UV detection with RPWELL

Objective

- Characterize single-stage and doublestage RPWELL-based configurations as UV detectors
- Achieve high gain and high detection efficiency under stable operation
 - Expose to events with higher primary charge (dynamic range)







UV detection with RPWELL

Results

• Polya pulse-height distribution observed with both configurations



A Polya Vs. Exponential

• Higher detection efficiency for equal average gain, noise and threshold



UV detection with RPWELL

Results

- Polya pulse-height distribution observed with both configurations
- Dual stage detector stable UV detection while exposed to 6 KeV x-rays background under a gain of $\sim 10^7$



Summary - The RPWELL

- Stable and efficient operation with single element (broad dynamic range)
 - Reduce the need for cascade configurations
 - Economic
 - Thin
 - Simple
- Upscaling is easy 50×50 cm² detector has been assembled and tested in the lab
 - Beam tests foreseen in July
 - Prototype sampling-element for (S)DHCAL foreseen end of the year
- Promising candidate for room-temperature UV imaging (e.g. in RICH)
 - Stable operation at room temperature with UV under 6 KeV x-rays background
- Under advanced study: Operation in Cryogenic temperatures:
 - Search for adequate resistive materials
 - Operation in Ne mixtures for reading light in noble-gas TPCs
 - Operation in vapor phase of noble gases

